

Update Description

The service pack CEREC SW 4.0.2 updates existing 4.0 and 4.0.1 installations and contains the following changes:

- **Preparation for CEREC Guide Data**
Adapted data handling in preparation for CEREC drilling guides
- **Bite Registration via Impression Material**
In the configuration menu (under Options), an option is available to add the image catalog type "Bite" to the "Add Catalog" field in ACQUISITION. Any pictures in this catalog are ignored if real images are taken in the catalog for the opposing jaw.
- **Language Option**
In the configuration menu (Options), the language of the user interface can be changed. The software must be restarted before changes take effect.
- **Material information Under "Change Block Size"**
In the optional step „Change Block Size“ within the MILL Phase, the block manufacturer and type are displayed.
- **All Restorations in Model Color**
If the analyzing tool "Color Restoration" is selected, inactive restorations are also shown in model color.
- **Contact colors on all Restorations**
Contacts are also shown on inactive restorations.

An existing CEREC SW 4.0 or 4.0.1 installation is required for this service pack. We strongly recommend this update to all 4.0 and 4.0.1 users.