CEREC[®] / inLab MC XL BULLETIN

12 S Bur Life Average

The 12 **S** Step Bur (P/N 6240167) and 12 **S** Cylinder Pointed Bur (P/N 6240159) were developed to provide the best balance of accuracy and bur life with the MC XL. Individual bur life on the MC XL can vary based on the size and type of block used, complexity of the milled restoration, amount of Dentatec used and water condition (cleanliness) and the milling mode (Fast or Normal).

Because of these variables, you can expect to get an average of 60 millings per package of Step or Cylinder Pointed burs. Please note that this is an average. There may be burs that need replacement after 4 to 6 millings while others will last 15 or more mills. This fluctuation is normal. However, if you experience consistent averages of less than 60 mills <u>per package of 6 burs</u>, please contact your Sirona authorized dealer for service.



Note: Life expectancy for the smaller 12 Step Bur (non-S) P/N 6260025 can be substantially lower than the 12 S Step Bur. Because of this, the use of this bur should be relegated to special cases where bur life would not be an issue and the smaller tip size would be of more benefit.





V1.0 Jun 2011